

# About the Tests

FoundationOne®Germline and FoundationOne®Germline More are next-generation sequencing (NGS)-based laboratory-developed tests designed to identify inherited genetic variants associated with hereditary cancer syndromes. These tests help guide risk assessment, clinical management, and potential targeted therapy selection for individuals with a personal and/or family history suggestive of hereditary cancer.

**FoundationOne Germline:** A 50-gene panel targeting clinically actionable hereditary cancer risk genes.

**FoundationOne Germline More:** A 154-gene expanded panel that includes clinically actionable and emerging cancer risk genes for patients requiring a broader risk assessment.

These tests are intended for use by qualified healthcare professionals to assess hereditary cancer predisposition based on clinical and family history. Results should be interpreted in conjunction with clinical findings.



## Analytical Performance

CHARACTERISTIC	ANALYSIS
Methodology	Hybrid-capture next generation sequencing (NGS)
Coverage	≥99% of targeted exonic regions at ≥50x depth
	Intronic coverage +/- 20 base pairs
	Deletion/Duplication detection as low as single exon resolution
Specificity and Sensitivity	Single Nucleotide Variants (SNVs): Sensitivity 99%, Specificity 99%
	Copy Number Variants (SNVs): Sensitivity >97%, Specificity >99%
	Insertions/Deletions (Indels): Sensitivity >97%, Specificity >99%
Reported Variants	Pathogenic, likely pathogenic, and variants of unknown significance (VUS)
Orthogonal Confirmation (if required)	Sanger sequencing and long-range PCR for selected variants, MLPA/qPCR for CNVs*
Reference Genome	GRCh37

\* Orthogonal confirmation includes MLPA (Multiplex Ligation-dependent Probe Amplification), qPCR (Quantitative Polymerase Chain Reaction) and long-range PCR

Both offerings also include analysis for the following regions/alterations of interest:

- *BRCA2* Portuguese founder mutation
- *MSH2* Boland inversion
- *PMS2* pseudogene resolution when applicable
- *CDKN2A*: p14ARF and p16 protein products
- Interrogation of promoter regions in select genes, such as: *GREM1*, *PTEN*, *TERT*, *TP53*



## FoundationOne Germline Gene List (50 genes)

APC	BRIP1	GREM1	MSH2	POLD1	SDHB	TSC2
ATM	CDH1	HOXB13	MSH3	POLE	SDHC	VHL
AXIN2	CDKN1B	MAX	MSH6	PTEN	SDHD	
BAP1	CDKN2A	MBD4	MUTYH	RAD51C	SMAD4	
BARD1	CHEK2	MEN1	NF1	RAD51D	STK11	
BMPR1A	EPCAM	MET	NTHL1	RET	TMEM127	
BRCA1	FH	MLH1	PALB2	RNF43	TP53	
BRCA2	FLCN	MLH3	PMS2	SDHA	TSC1	



## FoundationOne Germline More Gene List (154 genes)

ABRAXAS1	CDKN2A	EXO1	HRAS	PALB2	RET	TERT
AIP	CEBPA	EXT1	IKZF1	PALLD	RHBDF2	TGFBR1
AKT1	CEP57	EXT2	KIF1B	PAX5	RNF43	TINF2
ALK	CHEK2	EZH2	KIT	PDGFRA	RPS20	TMEM127
ANKRD26	CTC1	FAN1	LZTR1	PHOX2B	RUNX1	TP53
APC	CTNNA1	FANCA	MAX	PIK3CA	SAMD9L	TRIP13
ATM	CTRC	FANCB	MBD4	PMS2	SBDS	TSC1
ATR	CYLD	FANCC	MEN1	POLD1	SDHA	TSC2
AXIN2	DDB2	FANCD2	MET	POLE	SDHAF2	TYR
BAP1	DDX41	FANCE	MITF	POLH	SDHB	VHL
BARD1	DICER1	FANCF	MLH1	POT1	SDHC	WRAP53
BLM	DIS3L2	FANCG	MLH3	PRKAR1A	SDHD	WRN
BMPR1A	DKC1	FANCI	MRE11	PRSS1	SLC45A2	WT1
BRCA1	EGFR	FANCL	MSH2	PTCH1	SLX4	XPA
BRCA2	EGLN1	FANCM	MSH3	PTCH2	SMAD4	XPC
BRIP1	ELANE	FH	MSH6	PTEN	SMARCA4	XRCC2
BUB1B	ENG	FLCN	MUTYH	RAD50	SMARCB1	
CASR	EPCAM	GALNT12	NBN	RAD51C	SMARCE1	
CDC73	ERCC1	GATA2	NF1	RAD51D	SPINK1	
CDH1	ERCC2	GEN1	NF2	RB1	SRP72	
CDK4	ERCC3	GPC3	NHP2	RECQL	STK11	
CDKN1B	ERCC4	GREM1	NOP10	RECQL4	SUFU	
CDKN1C	ERCC5	HOXB13	NTHL1	REST	TERC	



## Limitations

This test detects germline (inherited) variants but does not assess for somatic alterations. Variants outside targeted exonic regions and deep intronic variants (>20bp from exon) may not be detected. Copy-neutral structural rearrangements (e.g., translocations) are not detected. Some variants, such as those in pseudogene regions, may require additional orthogonal confirmation. A negative result does not exclude hereditary cancer predisposition, as there may be undetected variants in genes not included in the panel or detectable by the test.

**For more details, visit:** [www.foundationmedicine.com/test/germline-testing](http://www.foundationmedicine.com/test/germline-testing)



## Regulatory Disclaimer

FoundationOne®Germline and FoundationOne®Germline More are laboratory developed tests that were developed, and their performance characteristics determined by Fulgent Genetics. FoundationOne Germline and FoundationOne Germline More have not been cleared or approved by the U.S. Food and Drug Administration.

**For clinical and technical support, contact:** 1-877-204-4319 or [FMISupport@fulgentgenetics.com](mailto:FMISupport@fulgentgenetics.com)

FoundationOne®Germline and FoundationOne®Germline More are laboratory developed tests that were developed, and their performance characteristics determined by Fulgent Genetics. FoundationOne Germline and FoundationOne Germline More have not been cleared or approved by the U.S. Food and Drug Administration. Specimen specifications are provided by Fulgent and are subject to change.